



VICART **Med**

Solutions for medical care

DISPOSABLE EQUIPMENT COVERS



DISPOSABLE EQUIPMENT COVERS

GCMedica offers a range of sterile single use equipment covers to maintain instant sterile field for theatre equipment in operating room. Our products include camera cover, operating light handle cover, C-arm drape and camera drape. Easy to use design provides cost effective and time saving solutions.





Camera Cover

Designed with lens and easy applying foam fixing ring. Available with square lens and round lens with different sizes for different camera.



Product code	Picture	Description
GC9006A03		Round len 100mm, foam base ID 78mm, length 12cm
GC9006A04		Round len 85mm, foam base ID 64mm, length 10cm
GC9006A05		Round len 57mm, foam base ID 45mm, length 12cm
GC9006B03		Square len 85mm, foam base ID 75mm, length 12cm



Light Handle Cover

Available with rigid light handle cover and flexible light handle cover.

Provides an instant sterile interface between non-sterile surgical lamp and the sterile surgeon and nurse.

We supply reusable metal adapters to standardize your lights to apply sterile light handle covers. Our adapters cover most operating room lights in market.









Product code	Picture	Description
GC9005B01		Flexible, green color with double rings
GC9005B02		Flexible, green color with rebar
GC9005B03		Flexible, green color
GC9005E01		Flexible, purple color
GC9005C01		Blue color
GC9005A01		Sponge fixing ring type, blue color
GC9005F01		Rigid, green color



Metal Adapter For Light Handle Cover

Standardize your lights to apply sterile light handle covers and quickly set up sterile field.

Various adapters fit most operating room lights in market.

Product code	Picture	Description	Applicable operating lights
GC1968A01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968B01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968C01		Fit for blue color light handle cover (GC9005C01)	MAQUET and STERIGIS AS lights
GC1968L01		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights
GC1968L02		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights
GC1968K01		Fit for green and purple flexible light handle cover (GC9005B01/ GC9005B02/9GC005E01/ GC9005S01)	XLED lights



Endoscopic Camera Drape

Endoscopic camera drape creates an instant sterile field around endoscopic camera heads and cabling. By eliminating the need for between-case sterilisation, these sleeve options can save money on costly repairs and prolong the life of equipment.





Quick-change Endoscopic Camera Drape

Buckle securing design allows quickly and easy scope change during procedures.

Closed buckle with high definition lens, acting as sterile barrier as well, offers camera full protection and extend its service life.

85mm broad ring reduces risk of contamination to the sterile field.



Product code	Picture	Description
GC1956A01		Size 160mm x 2.3m, Ring aperture 85mm



High Definition Endoscopic Camera Drape

High definition lens end ensures the HD image quality of camera. Closed end design offers camera full protection from contamination, acting as a reliable sterile barrier.

Allows quick and easy scope change during procedures without compromising sterility.

85mm broad ring reduces risk of contamination to the sterile field.

Product code	Picture	Description
GC1956B01		Size 160mm x 2.3m, Ring aperture 85mm



Elastomeric End Endoscopic Camera Drape

Elastomeric end can quickly wrap scope, helping easily and efficiently building up sterilized cover over endoscope. Drape is fully contained within the ring and can be simply pulled back over camera and cable. 85mm broad ring reduces risk of contamination to the sterile field.



Product code	Picture	Description
GC1956C01		Size 160mm x 2.3m, Ring aperture 85mm



Ring Folding Endoscopic Camera Drape

The drape permits the correct connection and operation of a surgical endoscope. Allows quick and easy scope change during procedures without compromising sterility. Open tip with two tapes.

Product code	Picture	Description
GC1956W01		Size 150mm x 2.5m, Ring aperture 85mm



C-Arm Drape

C-Arm drape is used for isolating X-Ray from the contamination source during the surgical operations.

It quickly builds up a sterile field for patients and doctors during procedures.

Clip fixing design for easy application.



Product code

GC9024M01

Description

Consist of 2 scope / fluoroscopy drapes
and one drape with 6 clips



VICART Med

Solutions for medical care

CONTACT US



**Chemin de Grottes 8D
1700 Fribourg - Switzerland**



+41 79 469 02 86



artur.cota@vicart.ch



www.vicart.ch